

DigiTUBEs

Disposable digestion tubes

DigiTUBEs are disposable digestion tubes suitable for use with commercially available graphite digestion block systems. They are manufactured with ICP-MS Grade Polypropylene and specially designed with a flat bottom to provide more direct contact with the base of digestion block wells, resulting in a faster heat-up time and more efficient sample temperature response times.

Every lot of **DigiTUBE** digestion tubes receives a Certificate of Volume analysis confirming the product's compliance with Class A Specification according to ASTM E288 and ISO 384 following our ISO/IEC 17025 accredited procedure.

For clean applications, **DigiTUBEs** are also available in shrink-wrapped bundles to protect against dust infiltration during storage and handling.



Key Features

- ▶ Ultra-low leachable metal content to minimize contamination risk
- ▶ Optional RackLock feature enables one-handed closure of screw caps
- ▶ Class A graduations for easy and accurate dilution
- ▶ Disposable, eliminating the need to clean the vessel
- ▶ Leak-free screw cap available in variety of colors

DigiTUBEs can replace 5 vessels in the lab - digestion beaker, graduated cylinder, volumetric flask, autosampler tube, and storage container. Save time and lower costs with the use of **DigiTUBEs**.

Bulk DigiTUBEs

| Description | Quantity | Part Number |
|------------------------------------|----------|-------------|
| 15 ml with RackLock, blue caps | 1620 | 010-515-261 |
| 15 ml with RackLock, no caps | 1620 | 010-515-262 |
| 15 ml without RackLock, blue caps | 1620 | 010-515-263 |
| 15 ml without RackLock, no caps | 1620 | 010-515-264 |
| 15 ml with RackLock, blue caps | 540 | 010-515-361 |
| 15 ml without RackLock, no caps | 540 | 010-515-363 |
| 50 ml with RackLock, blue caps | 750 | 010-500-261 |
| 50 ml with RackLock, no caps | 750 | 010-500-262 |
| 50 ml without RackLock, blue caps | 750 | 010-500-263 |
| 50 ml without RackLock, no caps | 750 | 010-500-264 |
| 100 ml with RackLock, blue caps | 300 | 010-501-261 |
| 100 ml with RackLock, no caps | 300 | 010-501-262 |
| 100 ml without RackLock, blue caps | 300 | 010-501-263 |
| 100 ml without RackLock, no caps | 300 | 010-501-264 |

Bulk Caps

| | Description | Quantity | Part Number |
|--|---|----------|-------------|
| Screw Caps for 15 mL <i>DigiTUBEs</i> | Screw Cap, Blue for 15 mL <i>DigiTUBEs</i> | 540 | 010-515-060 |
| | Screw Cap, Clear for 15 mL <i>DigiTUBEs</i> | 540 | 010-515-172 |
| Screw Caps for 50 mL <i>DigiTUBEs</i> | Blue Caps for 50 mL <i>DigiTUBEs</i> | 250 | 010-500-060 |
| | Yellow Caps for 50 mL <i>DigiTUBEs</i> | 250 | 010-500-130 |
| | Orange Caps for 50 mL <i>DigiTUBEs</i> | 250 | 010-500-140 |
| | Red Caps for 50 mL <i>DigiTUBEs</i> | 250 | 010-500-150 |
| | Clear Caps for 50 mL <i>DigiTUBEs</i> | 250 | 010-500-172 |
| Screw Caps for 100 mL <i>DigiTUBEs</i> | Blue Caps for 100 mL <i>DigiTUBEs</i> | 100 | 010-501-060 |

HandyPAKs

HandyPAKs include *DigiTUBE*s (50 mL or 100 mL) shrink wrapped to minimize contamination during storage or transport. Caps are provided in a separate bag. The **HandyPAK Plus** option contains pre-capped *DigiTUBE*s in convenient storage racks that are also wrapped for contamination control.

| | Description | Quantity | Part Number |
|----------------------|--|---------------|-------------|
| HandyPAK | 50 mL <i>DigiTUBE</i> s with RackLock, blue caps | 30 x 25 tubes | 010-500-362 |
| | 50 mL <i>DigiTUBE</i> s without RackLock, blue caps | 30 x 25 tubes | 010-500-364 |
| | 100 mL <i>DigiTUBE</i> s with RackLock, blue caps | 9 x 34 tubes | 010-501-362 |
| | 100 mL <i>DigiTUBE</i> s without RackLock, blue caps | 9 X 34 tubes | 010-501-364 |
| HandyPAK Plus | 50 mL <i>DigiTUBE</i> s with RackLock, blue caps | 20 x 24 tubes | 010-500-461 |
| | 50 mL <i>DigiTUBE</i> s without RackLock, blue caps | 20 x 24 tubes | 010-500-463 |
| | 100 mL <i>DigiTUBE</i> s with RackLock, blue caps | 24 x 12 tubes | 010-501-462 |
| | 100 mL <i>DigiTUBE</i> s without RackLock, blue caps | 24 x 12 tubes | 010-501-464 |



**Watch Glasses &
Storage Racks**

| Description | Quantity | Part Number |
|---|----------|-------------|
| Disposable Polypropylene Watch Glasses, 50 mL | 1000 | 010-500-081 |
| Disposable Polypropylene Watch Glasses, 100 mL | 500 | 010-501-081 |
| Plastic Storage Racks for 50 mL <i>DigiTUBE</i> s | 5 | 010-500-026 |
| Foam Storage Racks for 100 mL <i>DigiTUBE</i> s | 5 | 010-501-028 |